




Ekso Robotic Exoskeleton Clinical Research Summary of Findings

October 2017

Leslie VanHiel, PT, DScPT



Ekso Robotic Exoskeleton Research Update – 2017

Updated Numbers

Research with Ekso (referring to both first (Ekso1.1) and second (EksoGT) generation robotic exoskeletons) has grown in 2017. Currently, 62 studies (up from 55) have been completed or are ongoing internationally: 35 are completed, 27 are ongoing; and 6 are multi-center. These studies include the following diagnoses: spinal cord injury (SCI) (37), stroke (13), multiple sclerosis (3), traumatic brain injury (2), cerebral palsy (1), and multiple diagnoses (6). Total target enrollment is over 1,250 participants, with over 550 completed to date. Nine SCI studies and one stroke study have been published with several more in these and other diagnoses expected in the remainder of 2017. Ekso has been included in five published review articles¹⁻⁵ and one book chapter.⁶

Ekso and Walking

Ekso continues to demonstrate safety, progression of walking, and feasibility in a variety of diagnoses and settings (n=219).⁷⁻¹³ An additional study has supported users' learning curve of 15 sessions for walking with contact guard assist or close supervision (n=11).^{9,14} Preliminary data from an ongoing trial with Ekso in SCI¹⁰ and a systematic review¹⁵ supported previous findings¹⁶ that actual and perceived cardiorespiratory effort could be minimal or taxing during robot assist gait training walking, depending on level of volitional activity and walking speed. Thus, users can reap the benefits of walking with reduced cardiovascular fatigue and risk or with increased cardiovascular challenge.

Recently presented data further supported increased gait speed (n=31)^{9,10,14,17-20} and timed walking distance (n=23),^{9,10,17-19,21} as well as improved balance (n=18)^{19,22,23} while walking in Ekso, mostly for those with *complete* SCI. Of note, McIntosh¹⁰ showed 3 of 4 participants walked in the device faster than the threshold for limited community walking (0.44 m/s) for SCI. Ramanujam et al.^{20,23} determined that increased step and stride lengths and decreased stance time were highly correlated to increased gait speed in exoskeletons (n=11). Notably, Ekso was able to activate trunk and leg musculature below the level of injury in those with complete SCI (n=8)^{23,24}, adding to their postural control of weight shifts in Ekso. In those with *incomplete* SCI, Ekso targets muscle activation of proximal and posterior lower extremity muscles (n=14).^{23,25} Gait speed in this group measured outside of the device pre-post training has also been shown to significantly increase (n=29).²⁶⁻²⁸ Preliminary data²⁷ has shown gait speed improvement in participants with sub-acute to early chronic incomplete SCI (n=20), nearly meeting the minimal clinically important difference (MCID) for slow walkers (0.15 m/s)²⁹ and the threshold for limited community walking (0.44 m/s).³⁰ Significant improvements in timed walking distance (n=37)^{27,28,31} and walking balance (n=4)²⁸ have also been shown for this diagnosis.

Ekso has been shown to increase gait speed (n=30)³²⁻³³ and walking distance (n=23).³² in *stroke*. Molteni et al.³² described effects in people with sub-acute (n=12) and chronic (n=11) stroke after training 3x/wk for 4 weeks in Ekso. Both groups showed statistically significant improvement in gait speed and walking distance. Statistically significant improvements were also found for lower extremity strength and

ambulatory function in both groups, and for home/community ambulation in the sub-acute group.³² Finally, 5 out of 12 walked overground at baseline in the sub-acute group, whereas 9 walked at endpoint; 4 out of 11 walked overground at baseline in the chronic group, whereas 7 walked at endpoint.³² Preliminary data for acute inpatients with multiple diagnoses (n=6, mostly stroke), with mean days since diagnosis of only 16 days and half initially unable to complete the 10MWT, has shown average gait speed improvements after a mean of 7.3 ± 2.3 sessions with a mean of 416 ± 194 steps.¹¹ These improvements surpassed the MCID (0.16 m/s) and limited community walking threshold (0.40 m/s) for stroke³⁴ at 0.18 m/s and 0.44 m/s, respectively.¹¹ FIM transfer and locomotion scores also increased significantly.¹¹

Ekso and Balance

Human balance, or postural control, is the ability to “maintain, achieve, or restore [one’s] line of gravity within [one’s] base of support” to prevent a fall.³⁵ The balance measures BBS and TUG are highly correlated to walking in spinal cord injury (SCI) and stroke ($|r| = 0.81-0.93$).³⁶⁻³⁹

Several studies have shown that training in Ekso improves balance outcome measures in these diagnoses. A total of 11 participants with incomplete SCI, as well as those with complete SCI, decreased their time on the TUG while in the device, suggesting improved balance while supported by the exoskeleton.^{17,19} Outside of the device, significant improvements in TUG in incomplete SCI (n=4)²⁸ and in trunk control in stroke³² have been reported. Although not statistically significant, Flaherty et al.¹¹ reported an average improvement of over 5 inches in the Functional Reach Test outside the device (n=6).

The authors of the PanEuro study noted that the device relies on proper weight shifting to take the next step, and that training weight shifts should in turn train balance. Ramanujam et al.^{20,23} tracked the center of mass of participants with SCI during walking in Ekso (n=3), showing closer and more symmetric stability around midline²⁰, as well as “improved control ...of weight shifts and balance” after both high and low dosage training (n=7, range 3-100 hours).^{20,23} These combined changes demonstrated “a more efficient weight transfer strategy” and “improved dynamic stability”, resulting in a significant improvement in gait speed in the device. Motion capture in 5 participants with stroke showed weight shifts to be more centered with better symmetry of steps, while EMG showed more normal, distinctive lower extremity muscle bursts from the user.¹³ Improved midline alignment and excursion during standing outside the device after 4 weeks of training suggests better static stability as well.⁴⁰ Further studies of control of weight shifting during standing and walking are recommended, incorporating EksoGT’s SmartAssist software with real-time standing balance visual and audio feedback. This more normal gait pattern, more centered static and dynamic postures, enhanced strength, and higher dosage may enhance functional outcomes capacity during rehabilitation.

While the improvements in these outcome measures show increased *capacity* to walk, they do not necessarily translate into walking *activity* outside of the clinic. Low self-efficacy, or poor confidence to perform or complete a task successfully, has been studied as a limiting factor for walking post injury. French et al.⁴¹ studied 59 participants at least 6 months post stroke, finding that walking capacity results

significantly predicted self-efficacy results (Walk 12, Activities-Specific Balance Confidence (ABC) Scale, and Functional Gait Assessment (FGA)) and that self-efficacy results predicted walking activity and participation results. However, walking capacity results did not significantly predict walking activity or participation results. They determined that self-efficacy “mediated the relationship of performance-based [factors] with activity and participation” and suggest including self-efficacy measures in research and rehabilitation. In a similar study, in 55 participants at least 3 months post stroke, Danks et al.⁴² determined that a person’s walking capacity only explained 36% of the variability in his/her walking activity. A significant portion of the remaining variability was explained by 1) self-efficacy during balance tasks (ABC Scale) and during walking (Walk 12), and 2) the interaction of the ABC with the FGA. The statistical tests showed that balance self-efficacy was a significant predictive factor of walking activity and “moderated the relationship between walking capacity and walking activity”. The authors stressed that increasing self-efficacy through successes during rehabilitation may not only improve walking capacity, but also increase walking activity after discharge, thereby decreasing the secondary sequelae of immobility.

Ekso and Secondary Complications

Ekso has been reported to ameliorate a myriad of secondary complications following SCI and stroke, as well as increase quality of life. Two studies have shown significantly increased muscle volume and reduced intramuscular fat in chronic SCI (n=13).^{26,43} Lower extremity strength (n=28+)^{21,26,32,44} and bone density (n=6)^{43,44} have also been shown to improve after SCI and stroke. Participants have reported significantly reduced pain (n=18)^{14,18,19,44,45} and pain medication use (n=2)^{14,44} and reduced spasticity (n=32),^{8,14,21,45,46} after SCI, with some having trained 1-3x/wk for 5-20 sessions, but about half having walked in Ekso for only one session.⁴⁶ Participants with SCI have also reported improved bowel and bladder function (n=~18)^{14,21,44,46,47} and improved sleep (n=5).^{14,18} Preliminary data have so far shown 36 sessions of exoskeleton training (of which Ekso was one) reduced risk for cardiovascular disease via increased HDL cholesterol (n=~4), although LDL cholesterol, total cholesterol, nor triglycerides changed.⁴⁸ Notably, participants with SCI have reported significantly improved quality of life after 18 sessions (n=20).²⁷ Four participants with stroke reported they were more positive and motivated about their recovery when using Ekso.⁴⁹ Participants with SCI have reported improved mood and motivation (n=6)²¹, as well as reduced depression and improved body image (n=13)⁵⁰ with Ekso training. Users of Ekso with SCI have also described the psychological benefits of being upright and walking: appreciation of upright interactions with others, reduced anxiety and anger, and feelings of control, accomplishment, hope (n=8)⁵¹ and “a positive vision of the future” (n=13).⁵⁰

Ongoing Studies and Future Directions

The company-sponsored Walking Improvement for Spinal Cord Injury in Exoskeleton (WISE) randomized controlled trial began in 2016. Ten sites are included in this study to determine the effectiveness of a training protocol (3x/wk for 12 weeks) with EksoGT and SmartAssist software on gait, secondary complications, and PT and participant burden as compared to standard physical therapy for those with incomplete SCI. A similar RCT has begun in stroke this year with several participants already enrolled. Other ongoing studies include outcomes in walking, energy expenditure, and cognition in participants

with MS; bowel, bladder, muscle, and bone changes in SCI; blood flow changes in stroke and SCI; and psychological impact in multiple neurological diagnoses. A new, multi-center safety and feasibility study has just been launched to include deconditioned participants, whether in acute or long-term care rehabilitation. Using EksoGT and SmartAssist software, additional research is sought into changes in postural control and balance during walking, especially in those with midline alignment impairment and “pusher” syndrome after stroke.

Combining therapies often results in improved outcomes. Researchers at UCLA combined transcutaneous spinal cord stimulation (tSCS) with Ekso gait training for 5 sessions in one person with complete SCI to successfully evoke volitional lower extremity movement (outside of the device), perspiration, and cardiovascular control.⁵² Further study of adding tSCS and Busprione with Ekso training showed increased muscle activity and improved coordination patterns.⁵³ Brain stimulation before or during Ekso training via transcranial direct current (tDCS) or magnetic field (TMS) may be another way to augment outcomes. Despite the intervention, assessments of brain activity (EEG) and/or muscle activity/motor control (EMG) will be critical to future studies to show neuroplastic changes in the brain and spinal cord. We must understand these mechanisms underlying functional changes and relate them to participant characteristics of responders and non-responders to optimize EksoGT’s use in neurorehabilitation.

References

1. Louie DR, Eng JJ, Lam T, Spinal Cord Injury Research Evidence (SCIRE) research team. Gait speed using powered robotic exoskeletons after spinal cord injury: a systematic review and correlational study. *J Neuroengineering and Rehabil.* 2015;12:82.
2. Miller LE, Zimmermann AK, Herbert WG. Clinical effectiveness and safety of powered exoskeleton-assisted walking in patients with spinal cord injury: systematic review with meta-analysis. *Medical Devices: Evidence and Research.* 2016;9:455–466.
3. Nitschke J, Kuhn D, Fischer K, Röhl K. Reha-Technology comparison of the usability of the ReWalk, EksoGT and HAL exoskeletons in a clinical setting. *Orthopadie Technik.* 2014;22–26.
4. Contreras-Vidal JL, Bhagat NA, Brantley J, Cruz-Garza JG, He Y, Manley Q, Nakagome S, Nathan K, Tan SH, Zhu F, Pons JL. Powered exoskeletons for bipedal locomotion after spinal cord injury. *J Neural Eng.* 2016;13:031001.
5. Esquenazi A, Talati M, Jayaraman A. Powered exoskeletons for walking assistance in persons with central nervous system injuries, a narrative review. *PM&R.* 2017;9(1):46–62.
6. Jayaraman A, Burt S, Rymer WZ. Using Robotic Exoskeletons for Over-Ground Locomotor Training, a chapter in *Neurorehabilitation Technology*, Springer. 2016:493-511.
7. Nissen UV, Baunsgaard CB, Frotzler A, Brust AK, Ribeill C, Kalke Y-B, León N, G, Samuelsson K, Antepohl W, Holmström U, Marklund N, Glott T, Opheim A, Benito J, Murillo N, Nachtegaal J, Faber W, Gobets D, Biering-Sørensen F. Gait training in the EksoGT™ robotic-exoskeleton after spinal cord injury. Poster presentation. ISCoS Meeting; September 15, 2016.
8. Kolakowsky-Hayner SA, Crew J, Moran S, Shah A. Safety and feasibility of using the EksoGT™ bionic exoskeleton to aid ambulation after spinal cord injury. *J Spine.* 2013;S4:3.
9. Gorgey AS, Wade R, Sumrell R, Villadelgado, Khalil RE, Lavis T. Exoskeleton training may improve level of physical activity after spinal cord injury: a case series. *Top Spinal Cord Inj Rehabil.* 2017. Publication pending.
10. McIntosh, K. Exoskeleton assisted walking (EAW) in acute rehab following spinal cord injury (SCI): a safety and feasibility study. Platform presentation. ASCIP; September, 2017.
11. Flaherty V, Anderson C, Ball K, Furman K, Brady S. Use of the EksoGT Bionics exoskeleton robotic assisted mobility device with acute rehabilitation patients: preliminary results. Poster Presentation. APTA Combined Sections Meeting; February, 2017.
12. Russo A, Endersby K, Perret M, Kesten A, King M, Chervin K, Nolan K. A robotic exoskeleton to provide increased mass practice for gait training and its impact on discharge destination for individuals with acute stroke. Poster Presentation. International Stroke Conference; February 17-18, 2016.
13. Nolan K. Exoskeleton robotics in stroke rehabilitation. Platform Presentation. ACRM Meeting; November 2, 2016.
14. Kozlowski AJ, Bryce TN, Dijkers MP. Time and effort required by persons with spinal cord injury to learn to use a powered exoskeleton for assisted walking. *Topics in Spinal Cord Injury Rehabil.* 2015;21(2):110–21.

15. Lefeber N, Swinnen E, Kerckhofs E. The immediate effects of robot-assistance on energy consumption and cardiorespiratory load during walking compared to walking without robot-assistance: a systematic review. *Disabil Rehabil Assist Technol*. 2017 Oct;12(7):657-671.
16. Baunsgaard CB, Maher JL, Gerven JV, Palermo A, McMillian DW, Irwin RW, Biering-Sørensen F, Nash MS. Energy expenditure and cardiovascular drift effect during extended bionic walking. Poster presentation. ISCoS Meeting; September 16, 2016.
17. Stampacchia G, Rustici A, Bigazzi S, D'Avino C, Gerini A, Battini E, Franchini A, Tombini T, Mazzoleni S. Effects of exoskeleton gait training in SCI persons: pain, spasticity and endurance. Poster presentation. ISCoS Meeting; September 15, 2016.
18. Kressler J, Thomas CK, Field-Fote EC, Sanchez J, Widerström-Noga E, Cilien DC, Gant K, Ginney K, Gonzalez H, Martinez A, Anderson K, Nash MS. Understanding therapeutic benefits of overground bionic ambulation: exploratory case series in persons with chronic, complete spinal cord injury. *Arch Phys Med Rehabil*. 2014;95(10):1878–1887.
19. Sale P, Russo EF, Russo M, Masiero S, Piccione F, Calabrò RS, Filoni S. Effects on mobility training and de-adaptations in subjects with spinal cord injury due to a wearable robot: a preliminary report. *BMC Neurology*. 2016;16:12.
20. Ramanujam A, Spungen A, Asselin P, Augustine J, Canton S, Barrance P, Forrest G. Training response to longitudinal powered exoskeleton training for SCI. Abstract. International I Symposium of Wearable Robotics, 2016.
21. Umlauf S, Roth T, Schuurmans P, Müller F. Experiences with the exoskeleton EKSOGT BIONICS. Presented at DMGP Meeting; 2014.
22. Stampacchia G, Bigazzi S, Rustici A, Dini M, Corbianco S, D'Avino C, Gerini A, Battini E, Mazzoleni S. Psychological and physical effort during robotic exoskeleton assisted walking on treadmill and overground in SCI persons. Poster presentation. ISCoS Meeting; September 15, 2016.
23. Ramanujam A¹, Cirnigliaro CM², Garbarini E¹, Asselin P², Pilkar R^{1,3}, Forrest GF^{1,3}. Neuromechanical adaptations during a robotic powered exoskeleton assisted walking session. *J Spinal Cord Med*. 2017 April 20:1-11.
24. Alamro R, Chisholm AE, Lam T. Trunk muscle activation patterns during walking with robotic exoskeletons in people with high thoracic motor-complete SCI. Poster Presentation. ASNR Meeting; November 10-11, 2016.
25. Faber W, et al. The EksoGT skeleton - First clinical experiences with spinal cord injury patients and measurement of their energy expenditure in the Netherlands. Platform presentation. ISCoS Meeting; September 14, 2016.
26. Forrest G. Platform Presentation. ACRM Meeting; November 2, 2016.
27. Bonatti E, Baroncini I, Nalon C, Giovannini T, Falcone R, Mazolli N, Bonavita J. Locomotor training with exoskeleton EksoGT in patients with incomplete motor spinal cord injury in a hospital setting- preliminary results. Poster presentation. ISCoS Meeting; September 16, 2016.
28. Kern M, Chang S-H, Jennings J, Davis M, Francisco G, TIRR SCI Clinical Exoskeleton Group. Improving gait performance after SCI: a comparison between conventional PT and therapy using exoskeletons. Poster presentation. ASIA Conference, 2015.
29. Forrest GF, Hutchinson K, Lorenz DJ, Buehner JJ, VanHiel LR, Sisto SA, Basso DM. Are the 10 meter and 6 minute walk tests redundant in patients with spinal cord injury? *PLoS One*. 2014;9(5).

30. van Hedel HJA. Gait speed in relation to categories of functional ambulation after spinal cord injury. *Neurorehabil Neural Repair*. 2009;23(4):343-350.
31. Milia P, De Salvo F, Caserio M, Cope T, Weber P, Santella C, Fiorini S, Baldoni G, Bruschi R, Bigazzi B, Cencetti S, Da Campo M, Bigazzi P, Bigazzi M. Neurorehabilitation in paraplegic patients with an active powered exoskeleton (EksoGT). *Digit Med*. 2016;2:163-8.
32. Molteni F, Gasperini G, Gaffuri M, Colombo M, Giovanzana C, Lorenson C, Farina N, Cannaviello G, Scarano S, Proserpio D, Liberali D, Guanziroli E. Wearable robotic exoskeleton for over-ground gait training in sub-acute and chronic hemiparetic stroke patients: preliminary results. *Eur J Phys Rehabil Med*. 2017;Jan 24. E-pub ahead of publication.
33. Angacian G, Bialek A, Schwabe E, Cortes M, Rykman Peltz A, Climent A, Edwards D. Quantifying gait outcomes in chronic stroke using two robotic training protocols. Burke Rehabilitation and Research Student Presentation, 2015.
34. Perera S, Mody SH, Woodman RC, Studenski SA. Meaningful change and responsiveness in common physical performance measures in older adults. *J Am Geriatr Soc*. 2006;54(5):743-749.
35. Pollock AS, Durward BR, Rowe PJ. What is balance? *Clin Rehabil*. 2000;14:402-406.
36. Wirz M, Müller R, Bastiaenen C. Falls in persons with spinal cord injury: validity and reliability of the Berg Balance Scale. *Neurorehabil Neural Repair*. 2010;24(1):70-77.
37. Wang CH, Hsueh IP, Sheu CF, Yao G, Hsieh CL. Psychometric properties of 2 simplified 3-level balance scales used for patients with stroke. *Phys Ther*. 2004;84(5):430-438.
38. van Hedel HJ, Wirz M, Dietz V. Assessing walking ability in subjects with spinal cord injury: validity and reliability of 3 walking tests. *Arch Phys Med Rehabil*. 2005;86(2):190-196.
39. Flansbjerg UB, Holmbäck AM, Downham D, Patten C, Lexell J. Reliability of gait performance tests in men and women with hemiparesis after stroke. *J Rehabil Med*. 2005;37(2):75-82.
40. Michalek A, Janiszewska N, Cygon K, Maliszewski J. Reeducation of gait using exoskeleton EKSOGT GT in patient after ischemic stroke – case report. Poster Presentation. International Rehabilitation Days, February 12-13, 2015.
41. French MA, Moore MG, Pohlig R, Reisman D. Self-efficacy mediates the relationship between balance/walking performance, activity, and participation after stroke. *Top Stroke Rehabil*. 2016;23: 77-83.
42. Danks KA, Pohlig RT, Roos M, Wright TR, Reisman DS. Relationship between walking capacity, biopsychosocial factors, self-efficacy, and walking activity in persons poststroke. *JNPT*. 2016;40: 232-238.
43. Karelis A, Carvalho LP, Castillo MJ, Gagnon DH, Aubertin-Leheudre M. Effect on body composition and bone mineral density of walking with a robotic exoskeleton in adults with chronic spinal cord injury. *J Rehabil Med*. 2017;49:84–87.
44. Strausser K. Safety and efficacy of high-dosage use of exoskeleton in home environment for chronic SCI: a pilot study. Poster Presentation. ISCoS Meeting; September 16, 2016.
45. Stampacchia G, Rustici A, Bigazzi S, Gerini A, Tombini T, Mazzoleni S. Walking with a powered robotic exoskeleton: subjective experience, spasticity and pain in spinal cord injured persons. *Neurorehabil*. 2016;39(2):277–83.

46. Luard K, Martinelli L, Hobbs H, Faulkner J. A study exploring the clinical effects of a short-duration exoskeleton rehabilitation programme on key physiological markers in spinal cord injury. Poster presentation. ISCoS Meeting; September 16, 2016.
47. Hong EK. Patient-reported bladder management improvements after exoskeletal-assisted walking. Platform presented at ASCIP, Sept. 2017.
48. Knezevic S, Hong E, Asselin P, Kornfeld S, Cirnigliano CM, Gorman PH, Forrest G, Bauman WA, Spungen AM. Increased serum high density lipoprotein after 36 exoskeletal-assisted walking sessions. Poster presented at ASCIP, Sept. 2017.
49. Hammerton J, Wong S, Moki S, Kaduji A, Iftika S. What are the experiences of patients and clinicians from a spinal injuries unit using an EKSOGT® exo-skeleton device as part of an acute rehabilitation programme? Presented at Cybathon; October, 2016.
50. Milia P, De Salvo F, Caserio M, Cope T, Weber P, Santella C, Fiorini S, Baldoni G, Bruschi R, Bigazzi B, Cencetti S, Da Campo M, Bigazzi P, Bigazzi M. Neurorehabilitation in paraplegic patients with an active powered exoskeleton (EksoGT). *Digit Med*. 2016;2:163-8.
51. Stearns-Yoder KA, Dorsey Holliman BA, Penzenik ME, Hadidi DA, Lange SM, Hebert JR, Brenner LA. Exploring the psychosocial impact of EksoGT Bionics technology. Poster Presentation. ACRM Meeting; November 2-4, 2016.
52. Gad PN, Gerasimenko YP, Zdunowski S, Sayenko D, Haakana P, Turner A, ... Edgerton VR. Iron "ElectriRx" man: Overground stepping in an exoskeleton combined with noninvasive spinal cord stimulation after paralysis. *Proceedings of the Annual International Conference of the IEEE Engineering in Medicine and Biology Society, EMBS*, 2015, November:1124–1127.
53. Gad P, Gerasimenko Y, Zdunowski S, Turner A, Sayenko D, Lu DC, Edgerton VR. Weight bearing over-ground stepping in an exoskeleton with non-invasive spinal cord neuromodulation after motor complete paraplegia. *Front Neurosci*. 2017 Jun 8;11:333.

Other Ekso References

1. Molteni F, Gasperini G, Gaffuri M, Colombo M, Giovanzana C, Farina N, Lorenzon C, Scarano S, Cannaviello G, Guanziroli E. Neuromuscular pattern of the lower limbs of hemiparetic stroke patients during over-ground gait training: acute changes induced by a wearable exoskeleton. Poster presentation. Congregation of the Sisters Nurse Dell'addolorata.
2. Marrocco S, Aman S, Cornell S, Hsieh JTC, Fielding M, Wolfe DL. PRIME development: a data management and clinical decision support system to improve locomotor outcomes in persons with spinal cord injury. Poster Presentation. ISCoS Meeting; September 16, 2016.
3. Cornell S, Hsieh JTC, Fielding M, Marrocco S, Wolfe DL. Optimizing mobility outcomes across locomotor training modalities: clinical reflection during development of the PRIME algorithm – a case series. Poster Presentation. ISCoS Meeting; September 16, 2016.

